



# RISH SURGICAL

+92 3365240977

+92 3365240977

rishsurgicalinstruments@gmail.com

## Appendectomy Instrument Set



## **Backhaus Towel Forcep 13 cm**



The Backhaus towel clamp is a perforating clamp. It is used for grasping the tissue, securing towels or drapes, and holding or reducing small bone fractures. Backhaus Towel Clamp is mainly used to fasten drapes or towels to patients during medical procedures in order for them to stay on.

## **Metzenbaum Scissor 28 cm BL/BL, CVD**



Metzenbaum scissors or "Metz" are used for fine dissection and cutting. They are most frequently slightly curved.

## **Crile Forcep 14 cm 1:2 CVD and STR**



Crile Forcep 14 cm  
1:2 CVD



Crile Forcep 14 cm  
1:2 STR

Crile hemostatic forceps are atraumatic and non-toothed clamps used to clamp blood vessels, grab tissue or vessels. They are often known as a "snap."

## **Mayo Stille Scissor 15 cm BL/BL, CVD and STR**



Mayo Stille Scissor  
15 cm BL/BL, CVD



Mayo Stille Scissor  
15 cm BL/BL STR

Mayo-Stille Dissecting Scissor is a sturdy general surgical instrument designed for deep dissection of dense/tough tissues during general surgeries.

## **Mayo Hegar Needle Holder 16 cm**



Mayo Hegar Needle Holder is used to drive curved needles during surgeries. It has a broader jaw with a rounded tip. These needle holders have straight and cross-serrated inserts which provide greater precision during the suturing process. These needle holders are available in different sizes and variations.

## **Czerny Tenaculum Forcep 20.5 cm**



Czerny Tenaculum Forceps is a specialized surgical instrument that OB/GYN surgeons commonly use for grasping, holding and mobilizing the heavy internal female organs, in order to perform hysterectomies, oophorectomies, cervical tear repair and other procedures.

### **Schroeder Tenaculum Forcep 25 cm**



Schroeder Tenaculum Forceps are commonly used in surgery to hold blood vessels or steady the cervix and uterus. The forceps feature sharply pointed hooks at the tips of ring handles for comfortable and controlled use.

### **Lahey Seizing Forcep 20 cm / 8”**



Lahey Goitre Grasping Forceps assists surgeons to separate and dissect perithyroid tissues in goitre surgeries. For instance, they hold and cut the tissues related to the superior thyroid pedicle.

## Heywood Smith Biopsy Forcep 20 cm



Heywood Smith Biopsy Forcep is a specialized surgical tool that surgeons can use to grasp, hold and mobilize cotton wool swabs and other surgical dressings, in order to clean and prepare the operating site.

## Wertheim Cullen Forcep 21.5 cm



Wertheim Cullen Forcep is a forceps-like device that gynaecology and obstetrics specialists commonly use to manipulate heavier tissues, such as ligaments and other adnexa, during gynaecological surgery.

## Standard Dressing Forcep 13 cm, 25 cm



Standard Dressing  
Forcep 13 cm



Standard Dressing  
Forcep 25 cm

Standard Dressing Forceps is a hand-held device that surgeons in multiple areas can use to grasp and manipulate different types of wound dressings, in

order to perform wound debridement and remove necrotic tissue from the wound's surface.

### **Standard Forceps 1X2 T. 14.5 cm**



Tissue Forceps are used for in a various applications for holding and manipulating various types of tissues. Tissue forceps are commonly toothed in order to grip the tissue more firmly than serrated tips alone.

### **Scalpel Handle Standard Fig. 3 and No. 4**



Scalpel Handle  
Standard Fig. 3



Scalpel Handle  
No. 4

Intended use: To hold compatible blade for surgical incisions.

### **Yankauer Suction Tube 27 cm**



This tool is used to suction oropharyngeal secretions in order to prevent aspiration. A Yankauer can also be used to clear operative sites during surgical procedures and its suctioned volume counted as blood loss during surgery.

### **Middeldorpf Retractor 22 cm 20x22 mm**



Middeldorpf Retractor is a versatile surgical tool that allows surgeons to mobilize, retract and set tissues aside during plastic and maxillofacial surgeries. Usually, the tissues treated consist of delicate structures within the oral cavity, such as mucoperiosteal flaps, cheek, lips and the tongue.

### **Langenbeck Retractor 21 cm 30x11 mm**



Langenbeck retractors are commonly used within surgical procedures to separate the edges of an incision. They can also be used to hold back tissue or organs. We have both reusable, and single-use Langenbeck retractors available.

**Sterilization Container Stainless Steel 350x180x85 mm**

