


Embriotomy Instrument Set



Kelly Uterine Polypus Forcep 32 cm



Kelly Uterine Polypus Forceps is a specialized instrument commonly used for grasping delicate tissue within the abdominal cavity. Frequently used to grasp and remove polyps and fibromas from the used.

Smellie Perforator 27 cm




Smellie Perforator used during the procedure of Fetal skull decompression and misoprostol for cervix dilation. When the child had expired and the mother's life in danger, it was imperative that the dead fetus be removed immediately. Its overall length is 270 mm.

Teale Tenaculum Forcep 23 cm



The Teale Tenaculum Forceps offers multiple benefits for OB/GYN surgeons. Its principal use is to retract the uterus, fallopian tubes and blood vessels for performing hysterectomies and other surgeries.

Gauss 26 cm, Scalp Flap Forceps



Gauss Scalp Flap Forceps, Curved are used to grasp and clamp tissue in the uterus. The locking ring handle along with wide blades with grooves gets a firm grip on the tissue. The forceps are available in straight and curved.

Braun Decapitating Hook 32.5 cm




The Braun Decapitating Hook is used to decapitate the fetus after fetal demise inside the womb to aid vaginal delivery in obstructed labor, antepartum hemorrhage, uterine atony and maternal exhaustion. The instrument consists of a bulbed hook that is attached to a rod and then to a handle with an overall length of 30 cm. The obstetrician inserts the Braun Decapitating Hook into the vagina to lodge around the fetus and decapitate it to aid delivery. All of our Braun Decapitating Hooks are made with premium surgical stainless steel and can all be sterilized for repeated use.

Braun Cranioclast 42 cm, 1 Jaw



Braun Cranioclast is a specialized instrument that gynecological and obstetrics surgeons use to grasp, mobilize, and crush the fetal head, in order to ease the extraction of other fetal parts during embryotomy procedures.

Nagele Obstetrical Forcep 35.5 cm



Nagele Obstetrical Forceps is a specialized device that surgeons use to assist vaginal delivery, especially when the contraction force of the uterus is inadequate. This instrument may also be used during shoulder dystocia, breech delivery, and prolonged labor.

Kielland Obstetrical 42 cm



Kielland Obstetrical Forceps is a specialized device that obstetric surgeons use to grasp, mobilize, and pull the baby's head during complicated vaginal delivery, in order to minimize maternal and fetal mortality rates.

Sterilization Container SS 400x200x130 mm

