



RISH SURGICAL

+92 3365240977


+92 3365240977

rishsurgicalinstruments@gmail.com

Hecting Set



Spencer Scissor 9 cm, Normal



Spencer Stitch Scissor is used for dissection and cutting of internal and external sutures during a variety of surgeries including general, gynecological and vascular surgeries.

Standard Scissors 14.5 cm, BL/BL, CVD



Used for general cutting of sutures, gauze, and other materials; sharp blades used to cut tougher tissue and frequently used for dissection; blunt blades used for cutting materials such as tubing or gauze.

Halsted Mosquito Forceps 12.5 cm, STR/CVD



Halsted Mosquito
Forceps 12.5 cm, STR

Halsted Mosquito
Forceps, CVD 12.5 cm

With their fully serrated jaw and 5-inch working length, the narrow-pointed Halsted Mosquito forceps are typically used to clamp off smaller vessels that control the flow of blood.

Gross Maier Forcep 25 cm, STR



The Gross-Maier Dressing Forceps offers a broad range of surgical benefits. Its main purpose is to help surgeons in cleaning the blood and other obscuring fluids from vital surgical areas.

Fine Scissor 9 cm, SH/SH, Flat, STR



Fine scissors are used to cut and separate the soft tissues.

Mathieu Needle Holder 17 cm



Mathieu needle holders are used to hold needles for stitching during surgical procedures.

Forceps Narrow Dress 14.5 cm and Standard Forceps 1X2 T. 14.5 cm



Forceps Narrow
Dress 14.5 cm

Standard Forceps
1X2 T. 14.5 cm

Forceps are used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds.

Sterilization Container SS 350x180x85 mm

