




RISH SURGICAL

+92 3365240977


+92 3365240977

rishsurgicalinstruments@gmail.com

Herniotomy instruments Set




Backhaus Towel Forcep 13 cm




The Backhaus towel clamp is a perforating clamp. It is used for grasping the tissue, securing towels or drapes, and holding or reducing small bone fractures. Backhaus Towel Clamp is mainly used to fasten drapes or towels to patients during medical procedures in order for them to stay on.

Kocher Nippon Forcep 14.5 cm 1:2 CVD




Kocher Nippon Forcep is a specialized device that surgeons in different areas can use to manipulate blood vessels, clamp their lumen and control the blood flow, in order to achieve adequate haemostasis.

Collin Tongue Holding Forcep 17 cm



Collin Tongue Holding Forceps is a specialized device that ENT and maxillofacial surgeons can use to grasp, retract and statically hold the tongue, in order to expose abnormal tissue growths or lesions to excise.

Allis Baby Tissue & Holding FCP 13 cm 3:4



Baby Allis Tissue Forceps are widely used to grasp the intestinal tissues securely during surgical procedures. These forceps are nonreactive, healthy, and exceptionally balanced. These instruments are made from high-grade German surgical stainless steel and can be reused after sterilization.

Crafoord Sellors Ligature Forcep 24.5 cm



Crafoord Sellors Ligature Forcep is a specialized device that surgeons use to clamp blood vessels, especially small arteries, during a variety of surgical scenarios. The instrument is also useful to grasp disposable materials and soft tissues.

Crile Forcep Delicate 14cm STR



Crile Forceps Delicate are commonly used for clamping blood vessels or tissue before cauterization or ligation.

Standard Scissors BL/SH STR. 14.5 cm



Used for general cutting of sutures, gauze, and other materials; sharp blades used to cut tougher tissue and frequently used for dissection; blunt blades used for cutting materials such as tubing or gauze.

Foerster Forceps 25 cm, STR



Foerster forceps are essential for gynecologic surgery. Featuring oval tips with serrations, the Foerster is used to hold sponges during the procedure. The ratchet handle helps to hold the sponges firmly.

Mayo Hegar Needle Holder 20 cm



Mayo Hegar Needle Holder is used to drive curved needles during surgeries. It has a broader jaw with a rounded tip. They can be used for both human and animal surgeries. These needle holders have straight and cross-serrated inserts which provide greater precision during the suturing process.

Rochester Ochsner Forcep 16 cm, CVD




Rochester-Ochsner Forceps are hemostats that are helpful for clamping off blood vessels to avoid bleeding during procedures such as orthopedic surgery.

Standard Scissors 14.5 cm BL/BL, CVD



Used for general cutting of sutures, gauze, and other materials; sharp blades used to cut tougher tissue and frequently used for dissection; blunt blades used for cutting materials such as tubing or gauze.

Reverdin Needle 19.5 cm



Reverdin Needle are surgical devices that assist in passing surgical ligatures in narrow anatomical spaces.

Butterfly Grooved Director 14 cm



Butterfly Probe / Grooved Director is a specialized multipurpose device that surgeons in different areas use to probe cavities and wounds, enlarge fascia and isolate the tongue.

Scalpel Handle No. 4



Scalpel Handle No 4 is a hand-held device that surgeons in multiple fields use to carry and hold scalpel blades.

Forceps Narrow Dress. 14.5 cm and Standard Forceps 1X2 T. 14.5 cm




Forceps Narrow Dress. 14.5 cm



Standard Forceps 1X2 T. 14.5 cm

Forceps are used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds.

Potts Smith Tissue Forcep 21 cm, 1:2



Potts-Smith Tissue Forceps is a specialized instrument used for manipulating delicate tissues. It is used to access deep and narrow cavities as well as handle sensitive tissues during ophthalmic, ENT and plastic surgeries.

Collin Duval Forceps 20 cm 14 mm



Collin-Duval Forceps is a specialized instrument often used for gripping bulky tissues. Its innovative tips minimize damage while ensuring a secure grip. Commonly used in general, gynecological, orthopedic procedures as well as other specialized surgeries.

Farabeuf Retractor 12 cm



Farabeuf Retractor is one of the most popular hand-held retractors across the world, and it allows surgeons to grasp, retract and statically hold multiple types of tissues during procedures in a broad range of surgical specialties.

Parker Langenbeck 21 cm (set of 2)



Parker-Langenbeck Retractor is a multipurpose surgical device that assists surgeons in pulling back soft tissues and organs, and also in broadening wound edges or incisions during general surgery.

Retractor Depressor 250x25 mm



Sterilization Container Stainless Steel 350x180x85 mm

