



RISH SURGICAL

+92 3365240977

+92 3365240977

rishsurgicalinstruments@gmail.com

Minor Surgery Set



Standard Scissors BL/SH STR 14.5 cm



Surgical scissors are required for cutting skin, undermining the subcutis and deeper fascial layers, cutting sutures, and removing wound dressings

Mayo Scissor 15 cm STR/CVD



Mayo Scissor
15 cm STR



Mayo Scissor
15 cm BL/BL, CVD

Mayo scissors come various lengths. Their construction is extremely durable, and Mayo scissors are used to cut connective tissues, fascia and sutures, primarily. The tips are semi-blunt, and they are available in curved or straight styles.

Metzenbaum Scissor 18 cm BL/BL CVD



Metzenbaum scissors or “Metz” are used for fine dissection and cutting. They are most frequently slightly curved.

Crile Wood Needle Holder 15/18 cm



Crile Wood Needle
Holder 15 cm



Crile Wood Needle
Holder 18 cm

Crile-Wood Needle Holder is used to hold and guide small to medium size needles and suture materials. The Crile-Wood Needle Holder has a gently tapered, blunt tip and is generally more delicate than the Mayo-Hegar but slightly less delicate than Webster, Halsey, or DeBakey needle holders.

Foerster Forcep Serrated 18 cm



Foerster Forcep Serrated is a specialized surgical device that surgeons in every area can use for grasping, holding and mobilizing surgical dressings, especially Raytec sponges, for soaking up excess fluids.

Backhaus Towel Forcep 9 cm



Backhaus Towel Clamp is a popular device that surgeons in several fields use to statically hold surgical drapes to the patient's skin, in order to define the margins of the operating site.

Halsted Mosquito Forcep 12 cm STR



Halsted-Mosquito forceps are versatile, ratcheted, finger ring forceps used for clamping vessels and controlling blood flow and other clamping tasks in general surgical procedures.

Crile Forcep 14 cm, 1:2 CVD



Crile forceps are atraumatic and non-toothed clamps used to clamp blood vessels, grab tissue or vessels.

Spencerwells Artery Forcep 18 cm STR



Spencer-Wells Artery Forceps is a specialized device that surgeons in different areas can use to manipulate blood vessels, clamp their lumen and control the blood flow, in order to achieve adequate haemostasis.

Kocher Forcep 16 cm, 1:2 STR



Kocher Artery Forceps is a scissor-like device that surgeons in multiple fields use to grasp, mobilize and clamp bleeding vessels during wound closure, specialized surgery and other procedures.

Allis Baby Forcep 15 cm 4:5



Allis-Baby Forceps is a surgical instrument used for grasping tissues and occluding blood vessels during a variety of surgeries. It is commonly used during general, cardiovascular, gynecological as well as ortho and neuro surgeries.

Scalpel Handle Standard Fig 3



Scalpel Handle No 3 is a specialized device that surgeons in several areas use to statically hold scalpel blades, in order to incise the skin with maximum precision.

Forceps Narrow Dress. 14.5 cm, Standard Narrow Forcep 14.5 cm, Standard Narrow Forcep 18 cm, Standard Forceps 1X2 T. 14.5 cm



Forceps Narrow Dress. 14.5 cm



Standard Narrow Forcep 14.5 cm



Standard Narrow Forcep 18 cm



Standard Forceps 1X2 T. 14.5 cm

Forceps are used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds

Standard Tissue Forcep 18 cm 1:2



Tissue forceps are used in surgical procedures for grasping tissue. Often, the tips have "teeth" to securely hold a tissue. Typically tissue forceps are designed to minimize damage to biological tissue. These forceps are used for delicate tissue handling and traction during a procedure.

Volkman Bone Curette 17 cm Fig 0000



Volkman Bone Curette is a specialized device that orthopaedic surgeons use to gouge out abnormal bone growths and other hard tissues, in order to perform certain surgical procedures.

Butterfly Grooved Director 14 cm




Butterfly Probe / Grooved Director is a specialized multipurpose device that surgeons in different areas use to probe cavities and wounds, enlarge fascia and isolate the tongue.

Langenbeck Retractor 21 cm 30x11 mm



Langenbeck Retractor is a multipurpose surgical device that allows surgeons to perform retractions on delicate tissues, especially those found within the abdominal cavity, such as the intestines

Senn Mueller Retract 16 cm, 7x8.5 mm, SH



Senn-Miller Retractor is a multipurpose surgical device that allows surgeons to pull back and hold tissues away from the operating site during general surgery.

Sterilization Container Stainless Steel 350x180x85 mm

