



**RISH
SURGICAL**

+92 3365240977


+92 3365240977

rishsurgicalinstruments@gmail.com

Tracheotomy Instrument Set




Crile Forcep 14 cm 1:2 CVD



Crile forceps are atraumatic and non-toothed clamps used to clamp blood vessels, grab tissue or vessels.

Backhaus Towel Forcep 13 cm



Backhaus Towel Clamp is a popular device that surgeons in several fields use to statically hold surgical drapes to the patient's skin, in order to define the margins of the operating site.

Standard Scissors BL/SH STR 14.5 cm



Surgical scissors are required for cutting skin, undermining the subcutis and deeper fascial layers, cutting sutures, and removing wound dressings

Halsted Mosquito Forceps CVD/STR 12.5 cm



Halsted Mosquito
Forceps CVD 12.5 cm



Halsted Mosquito
Forceps 12.5 cm
STR


Halsted-Mosquito hemostatic forceps are versatile, ratcheted, finger ring forceps used for clamping vessels and controlling blood flow and other clamping tasks in general surgical procedures.

Mayo Scissor 17 cm, BL/BL CVD



Mayo scissors come various lengths. Their construction is extremely durable, and Mayo scissors are used to cut connective tissues, fascia and sutures, primarily. The tips are semi-blunt, and they are available in curved or straight styles.

Allis Forceps 15 cm 5x6



Allis Forceps is a popular surgical instrument that surgeons commonly use to grasp, statically hold or mobilize tissues of several kinds, in order to allow dissection, expose the surgical field or accessing a cavity.

Trousseau Trachea Trocar, Dilator 14 cm



Trousseau Tracheal Dilator assists surgeons to make incisions in the anterior tracheal wall during tracheostomy procedures. After widening the incision adequately, a tracheal tube is inserted in the stoma. Surtex Instruments recognize the vitality and the substantial nature of this procedure, and that it may be needed urgently to relieve acute airway obstruction. Maintaining the dilation of the tracheal opening using a dedicated device eases the insertion of the tracheostomy tube in the correct direction.

Jackson Trachea Retractor 15 cm



Jackson Tracheal Retractors are widely used to retract the delicate tissues during tracheotomy procedures

Metzenbaum Scissors 15/18 cm TC CVD



Metzenbaum
Scissors 15 cm
TC CVD

Metzenbaum
Scissors 18 cm TC CV

Metzenbaum scissors or “Metz” are used for fine dissection and cutting. They are most frequently slightly curved.

Mayo Hegar Needle Holder TC 18 cm



Mayo Hegar Needle Holder is used to drive curved needles during surgeries. It has a broader jaw with a rounded tip. These needle holders have straight

and cross-serrated inserts which provide greater precision during the suturing process.

Foerster Forceps 25 cm STR



Foerster Forceps is a specialized surgical device that surgeons in every area can use for grasping, holding and mobilizing surgical dressings, especially Raytec sponges, for soaking up excess fluids.

Forceps Narrow Dress 14.5 cm



Dressing forceps are used when dressing wounds. They hold gauze and other dressings. They may also be used during wound debridement to remove infected or necrotic tissue or debris from the wound. They may also be used for suturing.

Standard Forceps 1X2 T. 14.5 cm



Forceps are used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds

Adson Forceps 12 cm Serrated / 1X2 Teeth 12 cm



Adson Forceps
12 cm Serrated



Adson Forceps
1X2 Teeth 12 cm

Adson Forceps are used to handle tissues and dressings during surgical procedures. This instrument comes with a toothed tip and serrated surface.

Scalpel Handle Standard Fig 3



Scalpel Handle No 3 is a specialized device that surgeons in several areas use to statically hold scalpel blades, in order to incise the skin with maximum precision.

Frazier Suction Tube 19 cm Angled



Frazier Suction Tube is a specialized device that surgeons in multiple areas use to aspire fluids, debris and blood from multiple cavities, in order to achieve a greater view of the operating field.

Yankauer ASP Tube 23 cm Complete




Yankauer Suction Tube is a specialized device that surgeons in different areas use to dislodge particles, remove blood and fluids, and suction debris from the operating site.

Trachea Trocar & Dilator 18 cm



Tracheal dilation is performed to enlarge the airway in cases of subglottic stenosis and tracheal stenosis. Patients with these conditions complain of shortness of breath and dilation is a minimally invasive way of enlarging the airway.

Senn Mueller Retractor 16 cm 7x8.5 mm, SH



Senn-Miller Retractor is a multipurpose surgical device that allows surgeons to pull back and hold tissues away from the operating site during general surgery.

Chevalier Jackson Trachea Tube (50 mm D. 5 mm Fig 1) / (55 mm D. 6 mm Fig 2) / (75 mm D. 10 mm Fig 6) / (80 mm D. 11 mm Fig 7)



Chevalier Jackson
Trachea Tube 50 mm
D. 5 mm Fig 1



Chevalier Jackson
Trachea Tube 55 mm
D. 6 mm Fig 2



Chevalier Jackson
Trachea Tube 75 mm
D. 10 mm Fig 6



Chevalier Jackson
Trachea Tube 80 mm
D. 11 mm Fig 7

Chevalier-Jackson Tracheal Tube is a curved metallic assembly that can be placed in the tracheotomy opening to maintain airway patency. Metallic tubes

can be attracted by MRI machines and this must be explained thoroughly to prevent life-threatening accidents.

Sterilization Container Stainless Steel 350x180x85 mm

